

LOT EH0926

2026-11-05

Instrument	Parameter	Low	Normal	High
EDAN H60 H66 H68 H69 H60S H66S H68S H69S	WBC $\times 10^9/L$	3.47 $\pm$ 0.60	8.41 $\pm$ 1.00	19.19 $\pm$ 2.50
	NEU# $\times 10^9/L$	1.64 $\pm$ 1.0	4.46 $\pm$ 3.80	11.88 $\pm$ 7.60
	NEU% %	47.0 $\pm$ 10.0	53.1 $\pm$ 9.0	62.0 $\pm$ 8.0
	LYM# $\times 10^9/L$	1.3 $\pm$ 1.0	2.5 $\pm$ 1.5	3.5 $\pm$ 2.0
	LYM %	39.0 $\pm$ 10.0	29.6 $\pm$ 9.0	18.4 $\pm$ 8.0
	MON# $\times 10^9/L$	0.27 $\pm$ 0.20	0.73 $\pm$ 0.50	1.83 $\pm$ 1.10
	MON %	8.0 $\pm$ 4.0	8.7 $\pm$ 5.0	9.5 $\pm$ 6.0
	EOS# $\times 10^9/L$	0.20 $\pm$ 0.20	0.71 $\pm$ 0.50	1.90 $\pm$ 1.30
	EOS %	6.0 $\pm$ 5.0	8.5 $\pm$ 6.0	9.9 $\pm$ 7.0
	BAS# $\times 10^9/L$	1.00 $\pm$ 1.0	1.00 $\pm$ 1.0	1.00 $\pm$ 1.0
	BAS %	1.00 $\pm$ 1.0	1.00 $\pm$ 1.0	1.00 $\pm$ 1.0
	RBC $\times 10^{12}/L$	2.36 $\pm$ 0.30	4.48 $\pm$ 0.40	5.34 $\pm$ 0.50
	HGB g/L	62 $\pm$ 6	131 $\pm$ 8	168 $\pm$ 8
	HCT %	17.9 $\pm$ 3.0	39.2 $\pm$ 4.0	49.5 $\pm$ 5.0
	MCV fL	75.9 $\pm$ 5.0	87.5 $\pm$ 5.0	92.8 $\pm$ 6.0
	MCH pg	26.3 $\pm$ 2.5	29.4 $\pm$ 2.5	31.4 $\pm$ 2.5
	MCHC g/L	347 $\pm$ 30	336 $\pm$ 30	339 $\pm$ 30
	RDW-CV %	17.6 $\pm$ 3.0	15.3 $\pm$ 3.0	14.8 $\pm$ 3.0
	RDW-SD fL	50 $\pm$ 10	49 $\pm$ 10	50 $\pm$ 12
	PLT $\times 10^9/L$	58 $\pm$ 30	268 $\pm$ 45	519 $\pm$ 70
	MPV fL	10.1 $\pm$ 3.0	9.3 $\pm$ 3.0	9.2 $\pm$ 3.0
	PDW fL	12.8 $\pm$ 4.0	11.3 $\pm$ 4.0	11.2 $\pm$ 4.0
	PCT %	0.06 $\pm$ 0.06	0.25 $\pm$ 0.20	0.48 $\pm$ 0.20
	PLCR %	26.3 $\pm$ 10.0	20.5 $\pm$ 10.0	19.5 $\pm$ 10.0
	PLCC $\times 10^9/L$	15.3 $\pm$ 15	55.0 $\pm$ 20	101.2 $\pm$ 30

【NOTE】

- 1) The controls should be stored in the refrigerator (2~8°C). After opening, it will keep stable for 14 days when it is stored airtight at 2~8°C.
- 2) Please equilibrate the controls to room temperature (15~30°C) before using it.
- 3) Controls must be well mixed before use. Please mix gently to avoid cells rupturing or generating bubbles.
- 4) After using, put the controls back into the refrigerator (2~8°C) to prevent contamination and evaporation.